


With this Swiss Soft Laser technology a solution was found for treating tinnitus (ringing in the ears) and similar inner ear disorders, such as partial deafness due to tinnitus, ear pressure, morbus Ménière (dizziness), acute hearing loss and hearing distortion.



Tinnitus

Ringling in the ears

Laser treatment, no technical knowledge necessary

The Soft Laser System consists of the MedicLaser (5mW output power and 650nm wavelength) and the TinniTool® optical fibre (fibreglass optical fibre in the ear). The device was conceived so that treatment could be carried out at home (without technical knowledge). The effect of the MedicLaser + TinniTool® corresponds to photosynthesis in nature and is completely safe and painless. The laser beam even penetrates the deeper subcutaneous layers (Skin penetration: ca. 1.6 cm) and works as curative biostimulation directly

on the metabolism in the ear. More than 3000 affected persons have already been successfully treated with this system. Many disorders arise through a deficiency of ATP (adenosine triphosphate), the cell energy necessary for life. ATP is the key substance in our cells. Light is the only medicine that can increase the production of ATP. This needed light energy is stored in the 5,000 billions of cells. Photons (light particles) drive these cell power plants (adenosinetriphosphate). If the cell loses energy that is not re-

placed at the same time, this gradually leads to damage and finally to the death of the cell. Damaged sensory cells in the ear can also regenerate more rapidly through the concentrated administration of light energy. The concentrated light of the MedicLaser (650 nm) is optimally suited for that. The principle of TinniTool® is to direct the needed energy from the MedicLaser concentratedly and precisely to the ailing area. The laser treatment underwent successful long-term testing in hospitals and in physician practices. It was submitted to extensive medical tests and proven scientifically by various medical studies. The safety and effect of the ADisMa MedicLaser and TinniTool® were attested to and certified medically. The Soft Laser System (Battery operation: 2 x 1.5V, LR03 AAA) can either be purchased or rented for a period of three months.

MedicLaser + TinniTool®: EUR 288.-
or rental for three months: EUR 144.-

You will get more informations at:
DisMark GmbH
Rellikonstrasse 7, CH-8124 Maur
Phone: ++41 43 366 06 66
Telefax: ++41 43 366 07 66
eMail: info@dismark.ch
Internet: www.tinnitool.com